

# ACCUJECT™ REFRA



## Your choice for phakic IOLs

Make implantation of phakic lenses simple and cost efficient. Enjoy open eye access and a large loading chamber. Benefit from an injector that allows to push or screw without compromise. When it needs to be simple and cost effective - Accuject refra.

## ACCUJECT REFRA (for phakic IOLs / single-use / sterile)

Item Code	Device Name	Incision size (mm)*			IOL type
		wound-assisted	into-the-wound	into-the-bag	
AR3000	3.0 Injector	2.4	2.7 - 2.8	3.1	1-piece lenses
AR2600	2.6 Injector	2.2	2.5	2.8	
AR2400	2.4 Injector	2.1	2.4	2.7	
AR2200	2.2 Injector	1.9-2.0	2.2	2.5	

\* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques