

CE
2460

EYECRYL™ PLUS

Monofocal Hydrophilic Acrylic
Foldable Intraocular Lens



biotech

VISIONCARE

EYECRYL™ PLUS

Monofocal Hydrophilic Acrylic
Foldable Intraocular Lens


FEATURES

- 360 ° square edge
- 5 ° angulation ensures better contact with posterior capsule
- UV blocking Material
- Medical implantable grade hydrophilic material
- Packed in medical grade blister pack

SPECIFICATIONS

MATERIAL	Hydrophilic Acrylic CL UV	
OPTIC TYPE	Single Piece, 360° Square Edge	
OPTIC SIZE	6.00 mm	
OVERALL SIZE	12.50 mm	
ANGULATION	5°	
ACD	4.96	
REFRACTIVE INDEX	1.462	
RECOMMENDED ULTRASOUND A-CONSTANT	SRK-T 118.00	
RECOMMENDED OPTICAL A-CONSTANTS	HOLLADAY 1 SF: 1.42	HOLLADAY 2 SF: 1.45
	HOLLADAY 2 ACD: 5.2	SRK - T 118.40
	SRK - II 118.50	Hill-RBF: 118.30
	Barrett: 1.57	Olsen ACD: 4.20
	HOFFER Q ACD: 5.19	HAIGIS a0:0.940, a1:0.400, a2:0.100
	DIOPTRER RANGE	-10.0 D to +30.0 D (with 0.5 D steps from +1.0D to +30.0D)
IMPLANTATION SITE	Capsular Bag	
STERILIZATION	Steam	
SHELF LIFE	4 years from date of manufacture	

Model: 600



Scan For Biotech Connect App