


Delivery System for Hydrophilic IOLs

Brand Name	Model Name	Injector Name	Injector Color	Cartridge Type	Incision Size
Eyecryl Plus	600	Bio Hydroject BHSET-150	Green	Butterfly	2.8-3.0 mm
Eyecryl Plus Clear	HSAS600				



Note: Use only those injectors, which are supplied along with the hydrophilic segment IOLs. Do not inter change the injector with that of phakic and hydrophobic segment IOLs.

Delivery System for Hydrophobic IOLs

Brand Name	Model Name	Injector Name	Injector Color	Cartridge Type	Incision Size
Eyecryl Plus Natural HD	ASHFY600	Bio Hydroject BHSET-150	Purple	Butterfly	2.8-3.0 mm
Eyecryl ACTV+3	DIHFY600				
Eyecryl TORIC	HFY				
Eyecryl ACTV TORIC	HFYD				



Note:

- 1) Use only those injectors, which are supplied along with the hydrophobic segment IOLs. Do not inter change the injector with that of phakic and hydrophilic segment IOLs.
- 2) Eyecryl ACTV+3: DIHFY600, Eyecryl TORIC: HFY and Eyecryl ACTVTORIC: HFYD IOLs will have White color label to identify the stock compatible with 2.8 mm cartridge
- 3) This delivery system will be in place only until the stock of 2.8 mm compatible lenses get exhausted.