

Bio-Blue 90 Plus

Acid Blue 90 Ophthalmic Solution (0.025%w/v) With Polyethylene Glycol An Adjunct To Vitreoretinal Surgery

DESCRIPTION

- BIO BLUE 90 PLUS is a sterile, pyrogen free, non-toxic, isotonic and ready to use solution of Acid Blue 90 Plus (a non-fluorescent triphenylmethane stain) with effective staining properties.

INDICATION

- BIO BLUE 90 PLUS is indicated to selectively stain and visualize the internal limiting membrane (ILM), facilitating removal of tissue and reducing the risk of retinal damage. Before injection the vitreous chamber should be completely filled with suitable ophthalmic irrigating solution. Polyethylene Glycol prevents diffusion of this stain into other parts and due to its viscosity, it helps in controlled injection.

FEATURES

- BIO BLUE 90 PLUS selectively stains the internal limiting membrane (ILM), but not the epiretinal membrane nor the retina, making it easier to visualize the membrane for removal.

SUPPLY

Supplied in a 0.5 ml pre-filled syringe for easy opening and use.