

Bio-Blue Duo

Acid Blue 90 And Trypan Blue Ophthalmic Solution With Polyethylene Glycol An Adjunct To Vitreoretinal Surgery

DESCRIPTION

- BIO-BLUE DUO is a sterile, non-toxic, pyrogen free ready to use solution of Acid Blue 90, Trypan Blue and PolyethyleneGlycol. The syringe comes in medical grade blister packaging to ensure safety and protection

INDICATION

- BIO-BLUE DUO is a single injection designed for staining both the membranes, ILM as well as ERM and PVR membrane, without compromising the staining effect.
- It is injected in the vitreous cavity to clearly differentiate between the stained and unstained part of membranes and helps in easy removal of selected membrane.
- The combination with PEG 4% is found as stable dye and can also be injected in the eyes filled with BSS.
- Also it penetrates deep and stains only the targeted membranes or tissues.

FEATURES

- Stains ILM, ERM as well as FIVR membranes effectively.

SUPPLY

Supplied in a 0.5 ml pre-filled syringe for easy opening and use.