

Bio Blue Plus

Trypan Blue 0.15% w/v solution for Ophthalmic Vitrectomy Procedures

DESCRIPTION

- BIO BLUE PLUS is sterile, pyrogen-free solution of Trypan Blue 0.15% in saline buffer. It has selective affinity for Epiretinal Membrane (ERM). Therefore, it is useful as an aid in visualization of Epiretinal Membrane, facilitating its removal during Ophthalmic Surgical Vitrectomy procedures.

PROPERTIES

- Used as an aid during Ophthalmic Vitrectomy Procedures.
- Facilitates visualization and removal of Epiretinal Membrane and reducing the risk of retinal damage.
- It helps to identify and removal of posterior hyaloids remnants in Tractional Diabetic Retinopathy.

SUPPLY

Supplied in 1.0 ml glass vial as well as pre-filled syringe for easy opening and use.