

Bio Octane

DESCRIPTION

- BIO OCTANE is a sterile non-pyrogenic optically clear ultra-purified 100% fluorinated perfluorocarbon liquid used as an aid for retinal surgery.

PROPERTIES

- BIO OCTANE helps in the surgical management of retinal detachment
- Does not interfere with visualization of the retina
- BIO OCTANE has a high vapor pressure, which facilitates the removal of residual material remaining after aspiration

PRIMARY INDICATIONS

- Retinal detachment with giant tears.
- Retinal detachment with proliferative vitreoretinopathy.
- Traumatic retinal Detachment.
- Subluxated lenses.
- Retinal bleeding.

SUPPLY

Bio Octane is supplied in a glass vial as well as pre-filled syringe delivering 5.0 ml / 7.0 ml of sterile Perfluoro-n-Octane liquid. The product is packed in medical grade blister packaging and finally supplied along with information for use (IFU).