

Biosil / Biosil - F

A reliable and proven treatment for complex Retinal Detachments

DESCRIPTION

- BIOSIL range of products helps in the surgical management of retinal detachment. BIOSIL range of products acts as a mechanical tamponade and stabilizes retention force on the reattached retina.

PROPERTIES

- Non-Pyrogenic, clear and colorless
- Highly purified

PRIMARY INDICATIONS

- Proliferative Vitreoretinopathy (PVR)
- Retinal Detachment caused by Trauma
- Retinal Detachment with Giant Tears
- Retina detachments in the case of proliferative diabetic retinopathy (PDR)

Properties	Biosil	Biosil - F
Viscosity at 25°C ± 2°C	Between 950 cSt to 1050 cSt	Between 4500 cSt to 5500 cSt
Specific Gravity	Between 0.967 to 0.975	Between 0.967 to 0.975
Refractive Index at 25°C ± 2°C	Between 1.4000 to 1.4055	Between 1.4000 to 1.4055
Packing	10.0 ml PFS & Vial	10.0 ml PFS & Vial