


	EyeJet® CTR preloaded		
Model (REF)	13 RIGHT	13A RIGHT	13B RIGHT
Size (open)	12.3 mm	13.0 mm	14.5 mm
Size (compressed)	10.0 mm	11.0 mm	12.0 mm
Zonular Damage	up to 4 hours (120°)		
Incision	2.2 mm		
Material	BLUE PMMA		
Direction of Implantation	RIGHT (clockwise)		
Features	<p>The bent eyelet is designed to reduce a possible puncturing (bulging) of the capsular bag when the ring is being rotated into the bag.</p> 		

#### Advantages


- expands and stabilizes the capsular bag
- distributes the pressure to all zonular fibers and prevents one-sided pressure to single zonular fibers
- reduces post-operative capsular bag shrinkage
- facilitates centration of the IOL
- facilitates safe cataract surgery in eyes with weak or partially absent zonules
- the capsular tension rings have eyelets at both ends to facilitate insertion

CE 0483

PAGE 2/2

To use this product safely, please read "Instruction For Use" carefully before use.

 MADE IN GERMANY



## CAPSULAR TENSION RINGS

14 | 14C | 14A

	EyeJet® CTR preloaded		
Model (REF)	14 RIGHT / 14 LEFT	14C RIGHT / 14C LEFT	14A RIGHT / 14A LEFT
Size (open)	12.3 mm	13.0 mm	14.5 mm
Size (compressed)	10.0 mm	11.0 mm	12.0 mm
Zonular Damage	up to 4 hours (120°)		
Incision	2.2 mm		
Material	CLEAR PMMA		
Direction of Implantation	RIGHT (clockwise) / LEFT (counter clockwise)		

### Advantages

- expands and stabilizes the capsular bag
- distributes the pressure to all zonular fibers and prevents one-sided pressure to single zonular fibers
- reduces post-operative capsular bag shrinkage
- facilitates centration of the IOL
- facilitates safe cataract surgery in eyes with weak or partially absent zonules
- the capsular tension rings have eyelets at both ends to facilitate insertion

CE 0483

PAGE 2/2

To use this product safely, please read "Instruction For Use" carefully before use.

 MADE IN GERMANY

**Morcher GmbH**  
Kapuzinerweg 12 | 70374 Stuttgart  
GERMANY

Phone: +49 (0) 711 / 95 320 - 0  
E-Mail: [info@morcher.com](mailto:info@morcher.com)  
Web: [www.morcher.com](http://www.morcher.com)