

# TECHNICAL SPECIFICATIONS

MODEL	3FLA6E	3FLA6ET
MATERIAL	Hydrophobic Acrylic Foldable material	
OPTIC TYPE	Single Piece, 360° Square Edge with Refractive Aspheric Optic	
OPTIC SIZE	6.00 mm	
OVERALL SIZE	13.00 mm	
ANGULATION	0°	
ACD	5.28	
REFRACTIVE INDEX	1.48	
RECOMMENDED ULTRASOUND A-CONSTANT	SRK-T 118.60	
RECOMMENDED OPTICAL A-CONSTANTS	SRK - T 119.00	SRK - II 119.40
	HOLLADAY 1 SF : 1.74	HOLLADAY 2 SF : 1.79
	HOFFER Q ACD: 5.52	Barrett 1.88
	HAIGIS a0:1.309, a1:0.400, a2:0.100	
DIOPTER RANGE	+7.0 D to +30.0 D (with 0.5 D step) & 0.00 D to +6.5 D and +30.5 D to +40.0 D will be available on customization on case to case basis	
CYLINDER RANGE (Toric Version)	1.0 D to 1.5 D (with 0.5 D step), 1.5 to 6.0 D (with 0.75 D step).	
	6.75 to 15.0 D (with 0.75 D step) customization as per manufacturing possibility	
IMPLANTATION SITE	Capsular Bag	
STERILIZATION	Irradiation	

To choose suitable 3FLA6ET model, please logon to

