

Delivery System for Phakic IOLs

Brand Name	Model Name	Injector Name	Injector Color	Cartridge Type	Incision Size
Eyecryl PHAKIC	PKC120NH	Optiject BES-28P	Yellow	Butterfly	2.8-3.0 mm
	PKC125NH				
	PKC130NH				
	PKC135NH				
Eyecryl PHAKIC TORIC	PC120T				
	PC125T				
	PC130T				
	PC135T				

Note: Use only those injectors, which are supplied along with the phakic segment IOLs. Do not inter change the injector with that of hydrophilic and hydrophobic segment IOLs.

